

How to set up FMR sweeps based on static/conventional magnetometry data

Using SI units here: $B \text{ [tesla]} = \mu_0(H+M)$

For purposes of setting up initial FMR sweeps:

We can roughly estimate things easily by taking:

- $\gamma/2\pi \sim 28 \text{ GHz/T}$ (pure spin gyromagnetic ratio)
- expressing $\mu_0 H$ and $\mu_0 M$ in units of tesla (1 T = 10 kOe); sometimes we drop the μ_0 prefactor but the units of H or M in tesla mean that it is implied
- M_s can be calculated from $M(H)$ loop from SQUID, knowing film thickness and area, and using conversion:
 $M_s \text{ [tesla]} = M_s \text{ [emu/cm}^3] / 796$

Remember that the Kittel resonance condition reverts to the Larmor precession formula in the absence of anisotropies:

$$f = \gamma \mu_0 / 2\pi * H$$

and the discussion below is about solving the equation of motion for spins in the presence of various anisotropies relative to the direction of applied field:

- Demag field $H_D = -N_d M$ where N_d is the demag factor ($N = 1$ for thin film is assumed below)
- In-plane anisotropy field H_K^{in} which is the preference for moment to point along one direction in the plane of the film. This is usually very small in materials like CoFeB and NiFe.
- Perpendicular Uniaxial anisotropy field $H_K^\perp = 2K_U / \mu_0 M_s$
 - In our convention here, $H_K^\perp > 0$ tends to point moment out of plane
 - if H_K^\perp is large enough, namely larger than the demagnetization field, then $M_{\text{EFF}} = M_s - H_K^\perp < 0$, indicating we have PMA film
 - If $H_K^\perp = 0$ this implies an easy plane anisotropy: the material wants the moment to lie in the plane, but this does not specify any particular direction in-plane (that is specified by H_K^{in})
 - One can also have a tilted anisotropy in which H_K^\perp is not enough to completely bring the magnetization out of plane.

Applied field Out of Plane (OOP) of film:

$$f_{\text{RES}} = \frac{\gamma \mu_0}{2\pi} (|H_{\text{RES}}| + H_K^\perp - M_s) = \frac{\gamma \mu_0}{2\pi} (|H_{\text{RES}}| - M_{\text{EFF}})$$

So, an $H_K^\perp > 0$ adds to the applied field, while the demag field $H_D = -M_s$ needs to be overcome in order to bring the moment out of plane and produce a strong resonance ($f_{\text{RES}} > 0$). Note that H_K and M_s cannot be separated in OOP FMR measurements, so in our examples below we express the combined M_{EFF} .

- M_{EFF} is directly measured in SQUID $M(H)$ data along hard-axis:
 - PMA film: hard axis is in-plane, so in-plane saturation field = $H_K^\perp - M_s = -M_{\text{EFF}}$
 - IMA film: hard axis is out-of-plane, so OOP saturation field = $M_s - H_K^\perp = M_{\text{EFF}}$

Applied Field In Plane (IP) of film:

$$f_{RES} = \frac{\gamma\mu_0}{2\pi} \sqrt{(M_{EFF} + H_K^{in} + |H_{RES}|)(H_K^{in} + |H_{RES}|)}$$

H_K^{in} : in-plane anisotropy field; i.e., anisotropy between directions X and Y in the plane of film. $H_K^{in} > 0$ means it is collinear with external field.

The first term in the square root is the field seen trying to get the moment out of plane, and again M_{EFF} appears since it is the total perpendicular anisotropy (competition between demag and uniaxial anisotropies).

The second term in the square root is for the in-plane direction that is perpendicular to the external field.

We recover the familiar textbook Kittel equations for OOP and IP (Kittel's Intro to Solid State Physics, 5th Ed.) by ignoring H_K^\perp and H_K^{IN} .

Translating these equations to PhaseFMR software:



PhaseFMR	here	comments
IP M_s	M_{EFF}	Software doesn't know about H_K^\perp so it is included in reported M_s
IP H_K	H_K^{IN}	
OOP M_s	M_s	You can manually enter an $H_K (=H_K^\perp)$ in the software; we choose to leave that $H_K=0$ so that reported M_s is M_{EFF}
OOP H_K	H_K^\perp	

IP Example for IMA CoFeB film (50nm thick)

From SQUID: OOP saturation field = $\mu_0 M_{\text{EFF}} = 1.32 \text{ T}$

From FMR: $\mu_0 M_{\text{EFF}} = 1.267 \text{ T}$ ($\gamma/2\pi = 29.3 \text{ GHz/T}$ taken from Shaw(2015))

($H_{K_{\text{in}}}$ may be due to magnet remanence?)





IP Example for PMA CoFeB film:

From SQUID: IP saturation field = $-\mu_0 M_{\text{EFF}} = 0.154 \text{ T}$

From FMR: $\mu_0 M_{\text{EFF}} = \mu_0 M_s$ (in PhaseFMR) = -0.15 T

Reminder: in PhaseFMR below, H_k is actually H_k^{in} because we're doing IP analysis.




OOP Example for PMA CoFeB film (same one as above):



From SQUID: IP saturation field = $-\mu_0 M_{\text{EFF}} = 0.154 \text{ T}$

From FMR: $\mu_0 M_{\text{EFF}} = -0.124 \text{ T}$

(Agreement not good, may need to remeasure: we need to resolve this before publishing as an app note)



IP Example for Permalloy reference sample NiFe 10nm / Ta 5nm (sent in CryoFMR User Kit):
 Measured in-plane in “hard” direction (Sharpie line is perp. to H_{DC}); fitted $H_K < 0$ validates this
 From SQUID: OOP saturation field = $\mu_0 M_{EFF} = 0.91$ T (intersect lines from 0-0.6 T and 1.1-1.5T)
 From FMR: $\mu_0 M_{EFF} = 0.90$ T



OOP Example for Permalloy reference sample

