

APPLICATION TRAINING

**Data Collector
Data Viewer
HighScore and Highscore Plus**

XRD

August 2012

Copyright © 2012 by PANalytical BV, The Netherlands.

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system or transmitted in any form by any means electronic, mechanical, photocopying or otherwise without first obtaining written permission of the copyright owner.

This edition published in 2012 by:

PANalytical B.V.

Lelyweg 1, 7602 EA Almelo

P.O. Box 13, 7600 AA Almelo

The Netherlands

Tel: +31 (0)546 534 444

Fax: +31 (0)546 534 598

info@panalytical.com

www.panalytical.com

Contents

1. Safety Information	5
1.1 General Safety Information	5
1.2 General Safety Information	5
1.3 Safe Use of the Instrument	6
1.4 Diffractometer Enclosure	6
2. Startup/Shutdown Procedures	8
2.1 Turning the Instrument On	8
2.2 Data Collector	8
3. Basic Sample Preparation	12
3.1 The 'Back loading' Technique	12
3.2 Round Holders	12
4. Data Collection: Scan for Powdered Materials	20
4.1 General Considerations	20
4.2 Data Collector	24
4.3 X'Pert Industry	33
5. Data Analysis	40
5.1 Using DataViewer	40
5.2 Using HighScore	42
6. Routine System Maintenance	48
6.1 Basic System Maintenance	48
7. Specialized Applications	51
7.1 X'Celerator and PIXcel Detector Usage	51
7.2 Sample Changer Programs: Data Collector	63
7.3 Non Ambient Scans	66
7.4 CHC Humidity Chamber Measurements	69
7.5 Stress	77
7.6 Texture	82
7.7 Thin Film Grazing Incidence Analysis	84
7.8 Reflectivity	86
7.9 Data Base Conversion and Installation	90

Foreword

The information contained in this manual is intended to provide you with an introduction to basic system operation and user maintenance. The main focus is on how to collect an XRD powder scan and do an initial interpretation.

For further and more detailed information "User's Guides", "System Reference Manual" and "Quick Start Guides" for most packages are available as well as comprehensive information supplied in the on-line help for each package.

The topics covered in this manual are directed toward users of PANalytical's *Empyrean*, *X'Pert Powder* and *CubiX³* XRD systems.

The use of these instruments using the *Data Collector* and *X'Pert Industry* software is covered.

At publication of this document the latest versions of the software are:

- Data Collector Version 4.0(a),
- X'Pert HighScore Plus and X'Pert HighScore Version 3.0(d),
- X'Pert Industry Version 2.2(c),
- Data Viewer Version 1.4

Practical comments regarding more specialized applications are briefly covered in the supplemental chapters.

For the most current and detailed information refer to the "User's Guide", "System Reference Manual" or individual software "Quick Start Guides".

1. Safety Information

1.1 General Safety Information

The Empyrean, X'Pert Powder and CubiX³ X-ray diffraction systems conform to strict regulations regarding radiation and electrical safety. Please refer to the "System User Guide" for detailed safety information.

Refer to the supplied PANalytical Safety Manual and safety section of the User's Guide.

Some notable safety precautions:

- X-ray tube and detectors windows are made of beryllium. Do not touch these windows, and do not dispose of X-ray tubes or detectors with regular trash. PANalytical will take back our X-ray tubes and detectors for proper disposal. Look in XRD tube box for shipping label and instructions.
- If the shutter will not open, it is most likely due to part of the safety circuit not being met – doors improperly closed or safety lights not illuminated. Never try to override the safety circuits. Call the PANalytical Call Center on 1-800-279-7297 for technical support in help diagnosing the problem.
- Do not operate the system if the radiation safe panes on the door (on Empyrean and X'Pert Powder system) become damaged, cracked or broken. Turn power off and call the PANalytical Call Center on 1-800-279-7297 and open a service call.
- Before running samples verify there are no hazardous materials inside the instrument.

1.2 General Safety Information

1.2.1 Xsafe

The user is protected from unexpected exposure to X-rays by a safety system called Xsafe. This system comprises an electronic board which permanently monitors various switches throughout the system (for example checking that the doors are closed) to ensure that it is safe for X-rays to be generated.

A yellow LED next to the tube housing's shutter assembly illuminates when the shutter is open. If the safety circuit is broken during operation, a relay on the shutter assembly will automatically close the shutter.

1.2.2 Msafe

There is a further safety system in Empyrean, this is called Msafe. It ensures that users of the system are protected from any harm from system movement by ensuring that the goniometer and other modules cannot move while the doors are open.

Safety Information

Continued

1.3 Safe Use of the Instrument

During normal operation of Emyrean with all covers in position, no safety risks are posed to the operator. If the Msafe reset button needs to be pressed the lower rear panel must be removed in order to press the button. The following shows the hazardous areas to be avoided when that cover is removed.




Figure 1. Hazardous Areas – Rear

1.4 Diffractometer Enclosure

The upper part of the Emyrean Diffractometer Cabinet, is constructed from steel and is shown with the doors open in Figure 2. Those parts of the enclosure in the direct X-ray beam are 13 mm thick. Access to the inner part of the enclosure is via two fully interlocking doors at the front of diffractometer cabinet. The windows in the doors are made of lead glass with an X-ray absorption equivalent to 1 mm of lead.

Both of the doors actuate tongue operated guard locks in order to prevent accidents. As a safety measure, the tube housing shutter cannot be opened if the doors are not securely closed. An additional safety precaution is that the enclosure doors cannot be opened if the tube housing shutter is open, when the system is moving, or being initialized. The construction conforms with the most stringent X-ray safety standards: the dose rate is less than 1 $\mu\text{Sv/h}$ at 10 cm distance from the outside surface of the enclosure. Two lamps are fitted inside the enclosure. They can be switched on and off with the "LIGHT ON" pushbutton on the control and display panel.

WARNING:

If either of the windows is broken:


- Switch the system off immediately.
- Remove the safety key.
- Contact panalytical service.



Figure 2. Empyrean Diffractometer Cabinet - Upper Part

There are two power supply sockets (230 V, max. 4 A) mounted on the inside of the right-hand wall of the diffractometer enclosure.

A box with four LEMO connectors (1) and an Msafe indicator lamp (2) can be found on the lower right-hand side of the enclosure (shown below). You use these LEMO connectors to electronically connect motorized sample stages to the system. Connecting motorized sample stages to the system is described in the procedure for exchanging sample stages given in the user guide. The Msafe lamp lights if Msafe is activated. Resetting Msafe after it has been activated is described in the system user guide.



Please refer to the Empyrean Users Guide for more details on X-Safe, MSafe, MPM electronics, and (LEMO box)

2. Startup/Shutdown Procedures

2.1 Turning the Instrument On

This is the general procedure to start all XRD instruments. For more detail refer to the instrument specific quick start guides supplied with the instrument.

2.1.1 Preparation

1. Turn on the system computer (if off). Wait until Windows finishes loading.
2. Make sure that the chiller is on, water is flowing and the chiller is operating at the correct settings. (Make a note of the current pressure and temperature on the chiller gauges so you know what numbers to look for.)
Pressure = _____ PSI Temperature = _____ degrees F.

2.1.2 Turn the Instrument On

1. Use the 'Power On' button to start the XRD. The instrument will go through an initialization procedure which takes about minute.
2. If there is no reading on the kV/Ma windows check that the key is set to the on position.


2.1.3 Turn the Instrument Off


1. Turn the HV key on the front of the instrument to off.
2. Push the "Off" or "Standby" button".
3. Switch off the main power switch on the electrical box on the wall near the instrument.
4. Turn off the chiller approximately 1 to 2 minutes after the system has been powered off.

2.2 Data Collector:

2.2.1 Startup

- At the computer select the *Data Collector* from the programs menu or double click on the Data Collector icon on the desktop.
- When the Data collector Program opens go to the menu item "Instrument/ Connect"
- Select an appropriate configuration






Fine Calibration box will appear only if fine calibration offsets (FCO) are not zero. Press **'Yes'** button to continue using FCO or **'No'** button to clear the FCO to zeros.

2.2.2 Shutdown

1. From the Instrument Settings Control Panel window in *Data Collector*, double-click on the menu item labeled **'X-ray'** in the left side control window of Data Collector.
2. Uncheck the **'Generator on'** box and press **'OK'** to turn the generator off.
3. Push the Off or Standby button on the instrument control panel.
4. The Key may also be turned off. This is usually only done if there are radiation safety concerns.



2.2.3 Compressed air supply information

Source: House/building air or air compressor, if not available.

Purpose: X-ray Shutter operation, enclosure door latch, automatic beam attenuator, sample transport, etc.

	Specifications	Remark
Input pressure of instrument:	0.4 MPa to 0.5 MPa	4 to 5 atm (_____ , at installation)
Air consumption:	1.5 standard liter/sample At 0.4 MPa	Top open and close the shutter and door locks.
Lubrication:	No lubrication of mechanism at all	DRY air only!
Connecting hose:	5 m supplied with system	Inner diameter 8 mm. Outer diameter 14 mm.
Typical air compressor: 9425 901 91004, Jun-Air air compressor, Model 6X, 110 V 60 Hz		

Startup/Shutdown Procedures

Continued

2.2.4 Optional system add-on devices

Anton Paar DHS (domed hot stage) Compressed air Specifications

Pressure: 0.2 to 0.4 MPa (2 to 4 bar) (_____ , at installation)

Flowrate: 120 l/min at 0.2 MPa (2 bar) (_____ , at installation)

Anton Paar CHC (Humidity system) Compressed air Specifications

Pressure: 0.6 to 0.8 MPa (6 to 8 bar) (_____ , at installation)

Flowrate: 120 l/min at 0.7 MPa (7 bar) (_____ , at installation)

Dry air source for CHC

Purpose: Dry and wet air are mixed to yield desired Relative Humidity (RH%)

Source: House/building dry air or compressed air or nitrogen cylinder

Pressure: 0.2 to 0.6 MPa (2 to 6 bars or 30 to 90 PSI)

Quality: Class 2 – 5, Nitrogen gas

Install air preparation kit (i.e., 9425 011 91000 (US) or 9430 500 20521 (SC Almelo) to remove water and dust.

Important Maintenance: Drain water periodically, oil change, replace filter(s) on air prep kit.

2.2.5 X'Pert Industry:


Startup

On the Windows start menu or desktop icon open the *Industry* program.

The program will open and should communicate with the instrument



Manual hardware control can be obtained through the instrument button.



Followed by the '**Diffractometer Utilities**' button in the upper left corner of the Status window.

Select 'Generator' and input the desired power setting. Perform any necessary tube breeding procedures.



2.2.6 Shutdown

1. Use manual control to set the tube power down to 15KV and 5mA.
2. Push the Off or Standby button on the instrument control panel.
3. The Key may also be turned off. This is usually only done if there are radiation safety concerns.
4. Turn off the HV key, main circuit breaker, and chiller is necessary

2.2.7 Emergency Shutdown Procedure: All Instruments

If there is a need to immediately shut off the instrument and there is no time to power down the system using the software, press the "Off" or 'Standby' button on the front instrument console to immediately turn off the generator.

Note: This method is only recommended for emergency situations, as it can be stressful for the tube to be turned off abruptly.

2.2.8 EMERGENCY STOP Button/Switch (For Empeyrean and some X'Pert Powder (option) systems)

There is an EMERGENCY STOP button/switch on the front of the lower part of the cabinet. Should an emergency situation occur, press this button to immediately switch off the mains to the system. This switches the system off but leaves control and display panel and the safety circuits live. Before you can switch the system on again as described previously you must first release the EMERGENCY STOP button/switch by turning it clockwise until it springs out again.

3. Basic Sample Preparation

Diffraction measurements on powdered polycrystalline materials are the most common application. Common sample preparation techniques are presented here as an introduction. More comprehensive and detailed information is provided in the System's User's Guide.

3.1 The 'Back loading' Technique

Sample preparation procedures are critical in being able to obtain accurate and reproducible XRD results. Care should be exercised in order to avoid introducing errors resulting from factors such as:

- sample height displacement
- non-uniformity of the sample surface
- non-representative sampling
- contamination
- material loss
- alteration of composition due to over grinding, hydration, dehydration, or oxidation

3.2 Round Holders



Figure 3. Two-piece round sample holder



Figure 4. 9430 018 11271 Sample Holder and 9430 017 70101 Sample Preparation Kit.

Use the following guidelines to prepare a back-mounted sample:



Figure 5. Invert the holder ring and clamp it onto the sample preparation table. Spread the powder into the cavity using a spatula but do not pack or compress.

Basic Sample Preparation

Continued



Figure 6. Press powder with the Aluminum block. Remove excess powder with a knife blade or other straight edge.




Figure 7. Clean mating surfaces with small brush or edge of your thumb. Affix back plate.




Figure 8. Invert the prep table with holder attached while holding sample ring stationary, depress the release clamp by pushing on the button on the side and remove the sample holder. Check to ensure that the sample is level with the edge of the face of the sample holder.

3.2.1 Air sensitive sample holders



Figure 9. 9430 018 15401 Sample Holder Insert for Air Sensitive Samples, Zero Background Holder and Kapton Foil

Basic Sample Preparation

Continued




Figure 10. Shown here is a small amount of sample material onto the zero background sample holder evenly across the center.



Figure 11. A single sheet of X-ray transparent foil is placed over the zero background plate and the plastic clamping ring seals the foil onto the Sample Holder.




Figure 12. Cut away the excess foil outside the Sample Holder



Figure 13. Insert the Sample Holder into the PW1813/40 Sample Holder Ring




Figure 14. Sample Prepared in the 9430 018 15401 Air Sensitive Sample Holder and 9430 018 13401 Sample Holder Ring

Sample Preparation with Disposable Transmission Sample Holder Inserts



Figure 15. 9430 018 18401 Sample Holder Ring, Disposable Sample Holder Insert and Two Kapton Foils



Figure 16. Place the second Kapton foil over the remaining ring with the flat surface and clamp it with the ring that has the three raised bits



Figure 17. Clamping the Second Kapton Foil



Figure 18. Cut the excess foil extending beyond the circumference of the ring.

Basic Sample Preparation

Continued



Figure 19. Spread the Sample Material evenly across the center of the foil surface.



Figure 20. Clamp the Top Part onto the Bottom Part of the Insert



Figure 21. Place the insert into its position in the 9430 018 18401 Sample Holder Ring (a circular spring is not needed to clamp the insert into place).



3.2.2 Other Common Holders



Figure 22. 9430 018 17321 Si 'zero-background' plate. Scrape off excess material with a razor blade, leaving an even uniform surface of material level with the plate surface.



Figure 23. 9430 018 13XX1 Sample holder for Membrane filters, Metal plates, Pressed pellets.



Figure 24. 9430 018 12001 Steel Press Ring for pressed pellets in automated systems

4. Data Collection: Scan for Powdered Materials

4.1 General Considerations

4.1.1 XRD Scan Data Collection Strategy

It is a general practice to collect scans according to the following guidelines:

- Establish the relevant scan range including all major peaks. For example 5 to 100 degrees. The lowest angle in particular dictates the divergent beam slit settings and is described in more detail in the system user guide.
- For most inorganic materials 10 to 65 is suggested
- If clays are likely to be present 4 to 65 is suggested. Some specialized clay applications may require lower starting angles.
- For metals 30 to 120 (or the maximum angle) is suggested. If there are oxides or mineral type phases then 10 to the maximum angle is suggested.
- For organic materials 2 to 40 is suggested. For pharmaceuticals including inorganic component components 2 to 65 is suggested
- For organic materials transmission mode may be preferred due to preferred orientation or small sample amounts
- 5-10 data points above the full width of the half height of the maximum low angle peak is needed to accurately define a peak.
- This can be determined using the peak details tool in the DataViewer software.
- In order to obtain 1% relative counting statistical errors (for quantitative XRD analysis), collect approximately 10,000 total counts for the maximum low angle peak.

Many of optical configurations are possible for standard phase analysis/identification. The two main measurement types are focusing or Bragg Brentano optics and parallel beam measurements.

4.1.2 Focusing Measurements

Most common optics for a focusing measurement.

Application	Primary Beam Optics	Diffracted Beam Optics	Detector and Monochromator
Qualitative/quantitative analysis of bulk powders, flat solid samples, or small quantities on low background holder	Fixed or Programmable Divergence Slits	Fixed or Programmable Receiving and anti-scatter slits	<i>PIXcel</i> Detector <i>X'Celerator</i> Prop. Counter with Cu curved Monochromator

Typical configurations:

With *PIXcel* or *X'Celerator* Detector





Figure 25. Bragg-Brentano with *X'Celerator* or *PIXcel*

With Point Detector (Proportional)



Collection of a suitable measurement requires an appropriate selection of the following components:

Incident beam:

X-ray tube focus:	Line
Soller slit	Most common is 0.04 radians
Divergence slit	Depends on sample size, ¼ degrees is commonly used. Programmable slits may also be run in fixed area mode. PDS is used, fixed divergence slit settings of 1/32°, 1/16°, 1/8°, 1/4°, 1/2°, 1°, 2°, or 4° are possible.
Incident Anti-scatter slit	double the selection of the Divergent Slit. 1/2 degrees is commonly used
Beam Mask	Match to sample size, typically 15 mm
Beta-filter	Used with all detectors except when a monochromator is in place, typically not used on the incident side
Monochromator	Optional Germanium crystal for pure Alpha-1 geometry
Beam Knife	Optional accessory used for low angle measurements

Diffracted beam:

Receiving slit	With Point Detector, a common setting is 0.3 mm
Soller slit	Match with incident selection, typically 0.04 radians
Anti-scatter slit	Match to the selection of the Divergent Slit. For Linear Detectors refer to table below.
Beta-filter	Used with all detectors to remove beta radiation except when a monochromator is in place. Nickel (Ni) for a copper tube, iron (Fe) for a cobalt tube.
Monochromator	Used either to remove beta radiation or sample fluorescence. Used with point detector or linear detectors.
Detector	PIXcel, X'Celerator or Proportional Detector (Select one).



The following diagram is helpful in selecting the optimum divergence slit for measurements at the lowest scan angles.

These measurements conditions are independent of the software used to make the measurements.

Data Collection: Scan for Powdered Materials

Continued

Divergence Slit Selection Parameters




For a system with a fixed goniometer radius of 240 mm.
For other configurations, see <http://www.geol.uni-erlangen.de/html/xray/vcc.html> .

Table 1. Relationship between Empyrean System, Anti-scatter Slit Holder, Divergence Slit Size and Anti-scatter Slit Size for Symmetrical Scan Conditions

System	Empyrean	Empyrean with non-ambient Chamber
Anti-scatter Slit Holder	9430 034 94301	9430 034 94301 with 9430 034 94321
Divergence Slit Size (°)	Anti-scatter Slit Size (mm)	Anti-scatter Slit Size (mm)
4	13	No slit inserted
2	8.7	8.7
1	6.6	6.6
1/2	5.5	5.0
1/4	5.0	5.0
1/8	5.0	3.4
For alignment purposes	0.2	0.2

Note: The 3.4 mm slit is part of 9430 034 94321 Fixed Anti-scatter Slit. It is used in combination with 1/8° divergence slit and Anton Paar non-ambient chambers to measure at 2theta angles below 10°. More information is given in section ???????

4.1.3 Parallel Beam Measurements



Some possible optics for a parallel beam measurement.

Application	Primary Beam Optics	Diffracted Beam Optics	Detector and Monochromator
Thin film samples, bulk weak diffractors	Mirror	Parallel Plate Collimator	Prop. Detector, PIXcel or X'Celerator with 0D interface
Samples with rough or irregular surfaces, weaker diffracting thin films	Mirror	Parallel Plate Collimator	Prop. Detector, PIXcel or X'Celerator with 0D interface detector
Thin films or irregular samples containing Fe	Mirror	Parallel Plate Collimator	Prop. Detector, PIXcel or X'Celerator with 0D interface detector Flat graphite

Measurement Comments

Install the appropriate incident beam optic. For the mirror, insert the 1/2° slit. Use an appropriate beam mask based on the sample size with the Fixed or Programmable Divergent Slit or mirror. Insert the 0.04 radians Soller slits on the incident and diffracted beam optics.

Data Collection: Scan for Powdered Materials

Continued

4.2 Data Collector

4.2.1 Data Collector Basics

This section describes how to use Data Collector to set up a measurement. We presume that you have defined a user "Name" as described in Chapter 2 of the Data Collector Quick Start Guide [4].

4.2.2 Preparation

Switch the system on by pressing the "On" button on the control panel and then wait until the control panel shows: "30 kV" and "10 mA" indicating that the system is ready for use. If the power run up does not happen, refer to section at the beginning of this document.

4.2.3 Starting the software

Double-click on the Data Collector icon




Enter the user name and password: "My Name" and "password", and press . Default user name account and password are "User-1" and "Galaxy", respectively.



At this stage in the procedure go on-line (connect to the diffractometer)

- Select Instrument - Connect.



Select the appropriate configuration (in this example: "My Spinner").


If applicable - select the diffracted beam path with the Bragg-Brentano type optics and press .

A message window showing the 'assumed' status of the system is displayed:



In order to make sure that you obtain a good measurement, you must carefully check these assumptions. If these assumptions are correct, press and proceed with the next step. If they are not correct you must and then go to the tab(s) on the instrument control window containing the incorrect assumption and make the corrections.


Carefully look through the "Instrument Settings", "Incident Beam Optics" and "Diffracted Beam Optics" tabs to make sure that you have the correct parts mounted.



Data Collection: Scan for Powdered Materials

Continued

In the "Instrument Settings" tab double-click on the "Generator" path of the tree and enter "45 kV" and "40 mA":



Press and the system will power up to 45 kV and 40 mA.

4.2.4 Defining a Measurement Program

In this example we will first define a measurement program and then go on-line.

- Select File – New Program...




- Choose the type of program you want to define (in this example: "Absolute scan")



- and press

Defining the Measuring Program Parameters

You have just opened the “Prepare Absolute Scan” window:



- Select the appropriate configuration from the drop-down list. In this example we chose “My Spinner”. This will cause the following message to be displayed:



- Press
- In the “Configuration” frame, if applicable, select:

Diffracted Beam Path (these radio buttons only appear if your system has a double detector arm), choose the beam path with Bragg-Brentano focusing optics


- In the “Scan Axis” frame, select: Gonio (default)
- In the “Scan Properties” tab, select:
 - Continuous (default) Start angle (°) 20
 - End angle (°) 90
 - Step size (°) 0.020 (default)

NOTE: If you have a PIXcel or an X'Celerator detector in your system the step size is determined by the system according to the parameters of this detector (defined when you press the button).

Time per step (s) 1 The Scan speed and Total time are automatically calculated when you leave this field.

Data Collection: Scan for Powdered Materials

Continued



The next step is to define the hardware settings for the measurement that we have just defined.


Defining the Instrument Settings for the Measurement

Press **Settings...** to open a window in which we can specify the actual hardware settings that we are going to use.

To ensure repeatability and continuity please specify the components explicitly, (e.g. the description of the components as "Actual" should be replaced by the physical instrument settings).



Now select the stage and optics settings for this program, in this example we used the "Reflection-transmission spinner". Click on the Sample stage icon, select "Spinning" and if appropriate: the "Revolution time (s):" of "1".



In the “Incident beam path” branch of the tree we choose which of the items that we pre-selected that we want to use in this program, starting with the PreFIX module, in this example:

Progr. Div. Slit & Anti-scatter Slit:

Divergence slit:	Prog. Div. Slit (see Note following)
Usage:	Automatic
Irradiated length (mm):	10
Offset (mm):	0

Note: If you do not have a programmable divergence slit available to you, either select a fixed slit of 1°.

Mask:	Inc. Mask Fixed 10 mm
Soller Slit:	Soller 0.04 rad.

In the “Diffracted beam path” branch of the tree we choose which of the items that we pre-selected that we want to use in this program, starting with the PreFIX module (in this example: PRS/PASS):

Anti-scatter slit:	Prog. AS Slit
Usage:	Automatic
Irradiated length (mm):	10
Offset (mm):	0

Note: If you do not have a programmable anti-scatter slit available to you, select a fixed AS slit of 1°.

Receiving slit:	Prog. Rec. Slit (with a height of 0.1mm).
-----------------	---


Note: If you do not have a programmable receiving slit available to you, select a fixed receiving slit of 0.1/16degree.

Press to apply the settings to your program. Complete the creation of the measurement program by saving it: File – Save (in this example: “My Program”).

Check the “save” screen in the current version of the software (DC 4)

Data Collection: Scan for Powdered Materials

Continued




Press and close the Prepare Absolute Scan window by pressing .

4.2.5 Performing the Measurement

- Depending on your situation, mount the sample (Silicon sample in this case) as follows:

If you do not have a sample spinner in your system:	Mount the sample on the sample stage, close the enclosure doors and press "OK".
If you do not have a sample changer in your system:	Select the "Sample Stage" tab. Use the handle to lower the sample spinner platform, mount the sample, release the handle to bring the sample to the spinning position. Close the enclosure doors and press "OK".
If you have a sample changer, but it is positioned in the corner:	Select the "Sample Stage" tab. Uncheck "Lift Up", close the doors and press . Open the doors, mount the sample, close the enclosure doors, check "Lift Up" and press "OK".
If you have a sample changer, and it is ready to use:	Select the "Sample Changer" tab. Open the doors, put the sample into an empty magazine or monitor position. Close the doors, indicate to load the sample from the position you just loaded the sample into and press "OK".

Start the measurement program by selecting Measure – Program...




Highlight your program (in this example: "My Program") and press .

In the "File" frame enter:

Name: My Program_1.xrdml (default)
 Folder: C:\X'Pert Data\My Name (default)
 Comment: (in this example: "Example for Quick Start Guide").


In the "Sample" frame, give the sample a name:

ID: Standard Si Sample
 Name: Silicon Pellet
 Prepared by: select "My Full Name" from the drop-down list ().



Check that the enclosure doors are closed and press .

If any physical actions must be performed before the system can run the program, a list of these actions will be displayed.




• If there are any actions to be done, perform them and then press .


Data Collection: Scan for Powdered Materials

Continued

The scan starts and will take some time, depending on the program parameters (in this example: 1 hour). The scan is displayed as the measurement progresses:



Notice that the scale changes as the measurement proceeds.



This scan is automatically saved with the file name "My Program_1.xrdml".

You have now collected the data.


You can now use Data Viewer to view your results. A guide to using Data Viewer is given in the Explorer Add-ons Quick Start Guide.

Close Data Collector by selecting File - Exit and pressing

4.3 X'Pert Industry

4.3.1 X'Pert Industry

Collecting Basic Scans



Create a scan program. Go to 'Setup > Scan Program' and create a scan program such as the general program.






Figure 26. Input parameters for an X'Pert Industry scan

The diffractometer settings are found under "View Default Settings" and are assumed to have been set up correctly. By going into this option it is possible to check that the settings are appropriate for your applications. For example, if using programmable slits (PDS and PASS) with fixed parameters, make sure to change any values that may be input in the *automatic* PDS and PASS fields. Otherwise, the default automatic irradiated length setting (20mm) will override the fixed slit setting.

Data Collection: Scan for Powdered Materials


Continued




Select 'Save as' from the scan program window command line, name the scan program, and close the window.




Select '**Measure > Scan**' menu from the command line (not the 'Measure' button), Select the desired scan program name in the '**Scan program**' pull-down menu. Additionally, enter a '**File type**' (*.rd or *.xrxml), specify a '**Sample position**', and provide a '**File name**' for the scan. If a sample changer is in use samples can be loaded from the magazine(s) or the priority position (PP) underneath each magazine. Click on '**Start**' to begin the scan.



After the scan is complete, you can view and print graph from within *Industry* by clicking on '**Graphs > Scan**' on the toolbar button.




Browse and select a scan file, then press open display the scan.



Data Collection: Scan for Powdered Materials


Continued



To view the scan in *HighScore* for more detailed analysis open the program and from the command line, select '**File menu > Open**'. Load the file from the directory where it is saved. Make sure the file type specified on the browser page is "all files". Note that the default scan file location is "C:\PANalytical\X'Pert Industry\scans" but may be different on your computer.

4.3.2 Analysis Program


Most analysis in X'Pert Industry is made through an "Analysis Program"




Data Collection: Scan for Powdered Materials

Continued



Following is an example program which collects a scan over a Silicon Sample and returns the peak position and calculates the area of the silicon peak over a 0.8 degree range.



This screen selects the components to print out.



Analysis programs are run through the “Measure Sample” menu which opens a “Batch” where the measurements program, position on the sample changer and identity can be centered.



5. Data Analysis

5.1 Using DataViewer

5.1.1 Data Analysis and File Utilities

Data Viewer

Explorer Add-ons provides two additions to your Windows Explorer pop-up menu (View and Convert), info tips, thumbnail views of the contents of XRDML files, optionally, additional columns in the Windows Explorer view (described and shown in Chapter 2), and a data viewing facility with file conversion possibilities (Data Viewer).


Data Viewer provides facilities to enable you to view and work with measurement data. These extra functions are:

- Displaying a scan as graph(s)
- Displaying a measurement as an HTML report
- Displaying measurement data as ASCII text (valid for XRDML and other ASCII file types)
- Tool-like facilities
- Adding spectral lines to a graph.
- Adding peak parameters to a graph.

Data Viewer provides graphical functionality, reporting and conversion, integrated into the Windows Explorer environment. Search-match operations are supported through HighScore and highScore Plus.



For more information consult the "Explorer Add On Quick Start Guide". (Explorer Add-ons QSG.pdf)

Data Viewer can be started from a desktop icon or it can be executed automatically from Windows Explorer by double-clicking on most any datafile (e.g., *.xrdml, *.rd, or *.sd). Windows Explorer can be configured to display thumbnail views of all *.xrdml files in a given directory as shown below by selecting 'Thumbnails' from the View window.



Once Data Collector is started, files can be opened, viewed, overlain, and analyzed using a variety of tools that include, in part, size-strain and alpha-beta calculators, and tools for graph annotation, peak parameter determination, and report writing.


Most functions are selected through the menus or by using a mouse right click. Double click to unzoom.



Data Analysis

Continued

A file converter utility is also included in the Window Explorer add-on that is accessible by right-clicking on any *.xrdml file. Files can be converted to any of several formats, including PANalytical's older, binary (.rd) format and common ascii file types.



5.2 Using HighScore

5.2.1 Phase Identification with HighScore or HighScore Plus


HighScore is analysis software for powder diffraction measurements and phase identification. HighScore Plus adds Rietveld refinements and crystallographic analysis. It contains a lot of data visualization, data treatment and editing possibilities and reads almost all diffraction file formats. Parameter sets, user batches and a command line interface allow a high degree of customization and the automation of typical tasks. The user interface can be completely configured at will.

This discussion will cover a basic phase identification using an example file supplied with the software. More detailed information can be found in the quick start guide (XRD Quick Start Guides_ HighScore (Plus).pdf) and in the online help.

5.2.2 Routine Operation Using a Scan of a Mineral Mixture

Open *HighScore* and load the file '*minmix.rd*':


File > Open > C:\Program Files\PANalytical\X'Pert *HighScore*\Tutorial\ *Minmix.rd*'



Data Analysis


Continued

HighScore is extremely customizable and can have very different appearances. If it looks very different the default appearance can be set through the "View/Desktop/Phase-ID" menu. This can be useful if the Lists Pane is not visible.



To complete a phase identification

- Press the 'IdeMin' button in the lower left-hand corner of the screen to execute a batch search-match analysis of the data.
- Once executed, matched phases are identified and highlighted in the diffraction pattern, color-coded to a legend that identifies these *accepted* phases in the upper right-hand corner of the display. A larger *candidate* list is also displayed in the lower right-hand corner of the screen.
- The value of a 'Score' measures the strength of a phase match, the higher the score the better the match.
- If you see the message "Batch step not executed" disregard it as this sample has been collected in fixed divergence slit mode rather than automatic divergence slit mode and a conversion is not required.




The analysis functions embedded in the 'IdeMin' batch include background definition, peak search, and a search-match scheme from the ICDD minerals subfile.

Most screens and options can be modified using a right Mouse selection.

Data Analysis

Continued


A word document report can be created using the "Reports/ Create Word Report/Default" menu.



In addition to the power of rapid pattern processing, High Score can calculate the approximate weight fractions of each phase via a standard-less SemiQuant function. SemiQuant uses the RIR's (Reference Intensity Ratios) from the ICDD database in its calculation. If each matched phase has an RIR value then quantitative information is calculated.

No.	Visible	Ref. Code	Compound Name	Chemical Form.	Score	Scale	Display Co.	Quality	Subfiles	Crystal System	SemiQuant	Database ID
1	✓	WARR 88...	Quartz low	O2 Si1	99	8.998	Blue	C	ICSD ...	Hexagonal	28	D/Progra...
2	✓	WARR 88...	Hematite	Fe2 O3	55	8.412	Line	C	ICSD ...	Hexagonal	15	D/Progra...
3	✓	WARR 88...	Aluminum Oxide	Al2 O3	57	8.172	Gray	C	ICSD ...	Hexagonal	21	D/Progra...
4	✓	WARR 88...	Calcite	Ca1 Ca1 O3	48	8.618	Mbr	C	ICSD ...	Hexagonal	27	D/Progra...
5	✓	WARR 88...	Calcium Fluoride	Ca1 F2	46	8.428	Aqua	C	ICSD ...	Cubic	42	D/Progra...

Example of a SemiQuant analysis for a mineral mixture.



6. Routine System Maintenance

6.1 Basic System Maintenance

6.1.1 X-ray Tube Maintenance

The most current tube information is in the XRD Tube Instruction Manual included with every new tube.

Water Flow Monitoring

The X-ray tube in the diffractometer typically uses up to 1.8 kW for Copper and Cobalt or 3 kW for Molybdenum. Most of this power is converted to heat. The tube is cooled by means of an external air-cooled or water-cooled chiller during production of X-rays. An uninterrupted water flow of at least 3.7 liter (0.98 gallon) per minute with a maximum inlet temperature of 35°C (95°F) is required for effective anode cooling. A safety switch in the cooling system will switch the high tension generator off if the water flow falls below 3.7 liter (0.98 gallon) per minute.

Customer Service Engineers typically check, clean or replace the filters in the system during a Preventative Maintenance (PM) visit. Users should monitor the flow parameters of their systems between preventative maintenance visits as the filters can become clogged with debris resulting in reduced flow.

- Appropriate cooling water conditions are found in the Tube manual section 7.3.3 (X-ray Diffraction Empyrean Tube – Instruction Manual, June 2010, 4022 196 7141)
 - i) Flow: 3.7 – 6.0 liters (0.98 – 1.58 gallons) per minute.
 - ii) Max. Pressure: 0.8 MPa (8 bars or 116 PSI).

To monitor the water flow during instrument operation, use the following procedures:

Customizing the Status Bar in *Data Collector* to Monitor Water Flow

In *Data Collector*, connect to the system with a valid configuration. From the command line menu, go to '**Customize > Status Bar**' and select one of the five drop-down boxes to readout the water flow. The current water flow will now be displayed in the bottom of the *Data Collector* window.

Water Flow Monitoring in *Industry*

On the main toolbar, press the '**Instrument**' icon. Water flow conditions will be displayed in the lower right-hand corner of the instrument schematic.

Checking and Cleaning/Replacing Filters

Water filters are built into the cooling system to protect the system from dirt deposits. There are three filters in the water flow path:

- a **screen filter inside the chiller** that can be removed and gently cleaned with a brush,
- a **'whole-house' filter cartridge** between the chilled water supply and the main water inlet of the instrument. This filter should be replaced every 4 to 6 months. The User's Guide describes how to change this filter.
- a small **screen inside the water inlet of the X-ray tube housing** is located just before the water reaches the rear of the tube anode which can be removed and cleaned. The procedures for changing or cleaning the water filter in the X-ray tube are described in Chapter 5 of the XRD Tube Instruction Manual).

PreFIX Optics and Stages

When mounting PreFIX (Pre-aligned Fast Interchangeable X-ray) modules inspect and clean with supplied chamois and oil

Tube Breeding

During tube breeding (conditioning), pre-programmed voltage and current combinations are sent to the instrument in order to clear any impurities from the X-ray tube and to achieve a stable vacuum. Tube breeding can be performed by means of *Data Collector*. For more information on the tube breeding procedure refer to the XRD Tube Instruction Manual.

Please note that the XRD tube manual refers to "breeding" as "conditioning"


Likewise “Normal” is “Long” and “Fast” is “Short”

	Time since last switch OFF:	Time since last switch OFF:
	< = 100 hrs	>100 hours
Tube has NOT been removed from the system	Breeding is NOT required	Execute FAST (SHORT) breeding program
Tube has been removed from the system	Execute FAST (SHORT) breeding program	Execute NORMAL (LONG) breeding program

‘Normal’ Tube Breeding in *Data Collector*

This function requires that the user has diffractometer access level 2 privileges. To perform a NORMAL (LONG) tube breeding in *Data Collector*:

- Select from the menu ‘Tools > Tube Breeding, > Normal’.
- Alternatively, double-click the left mouse button anywhere in the ‘Instrument Settings’ tab of the Control Panel to open up the following window, and then click on ‘Breed’.



Instrument settings menu in *Data Collector*.

A normal tube breeding procedure typically takes about 15 minutes. Note that the ‘unusual’ kV/mA combinations used by this procedure are not harmful to the instrument. The program should be permitted to run to completion. The system will be unavailable during this process.

‘Fast’ Tube Breeding in *Data Collector*

To perform a ‘fast’ (SHORT) tube breeding in *Data Collector*, go to ‘Tools > Tube Breeding, > Fast’. It will take approximately 3 minutes to complete the procedure.

‘Normal’ Tube Breeding in *Industry*

To perform a ‘normal’ tube breeding in the *Industry* program, go to menu ‘Utilities > Diffractometer Utilities > TubeBreeding’. Again, allow the program to run to completion.

‘Fast’ Tube Breeding in *Industry*

To perform a FAST tube breeding in *Industry*, select ‘Utilities > Diffractometer Utilities > TubeBreedingFast’.

6.1.2 Standby mode

For maximum tube life and to eliminate the need for tube breeding during short (e.g., overnight or weekend) periods of no use, maintain the generator at 45kV and 20mA. The instrument will be ~immediately available for use after reconnecting to the software.

Routine System Maintenance

Continued

6.1.3 Water quality

Water quality is described in Appendix C of the XRD Tube Instruction Manual.

Table 2. Cooling Water Specifications

Parameter	Specification	Remarks
Water quality	Recirculating cooling system preferred, filled with drinking water.	
Water temperature	17 to 35 °C, and > dew point	See note below.
Water temperature fluctuation	< 2 °C/h sinus wave	
Water pressure	0.1 to 0.8 MPa	High pressure (at least 0.1 MPa overpressure behind the tube) is recommended.
Water flow	4 to 6 l/min	Stable and uninterrupted.
pH	7.5 to 9.5	Add Na ₂ CO ₃ to increase pH value.
Hardness	1 to 4 DH	DH = German Dedrees. Value can be controlled using a filter.
Conductivity	1 to 25 mS/m	To increase conductivity, add drinking water. To decrease conductivity, renew whole water volume with fresh drinking water.
Water filter (supplied with system)	50 µm filtering	Wall mounted.
Connecting water hose	10 m	
Connecting water pipe (not supplied with system)	Outer diameter 15 mm	Use the same materials as much as possible. Preferable syntjetic, otherwise stainless steel. To prevent algae non-transparent (opaque) materials should be used

The water pH should be 8.5 ± 1 and electrical conductivity should be 13 microSiemens/cm. Reference Emphyrean or PRO user guide for more information. Add note that chiller must be shut down (or waters valves closed to XRD system) must be when XRD system is power off to extend XRD tube life.

7. Specialized Applications

7.1 X'Celerator and PIXcel Detector Usage

7.1.1 X'Celerator Detector

The X'Celerator is a solid state linear X-ray detection system that can be used for powder diffraction applications. The system is based on Real Time Multiple Strip (RTMS) technology that permits significantly shorter measurement times than conventional detectors while maintaining or improving signal-to-noise characteristics and peak resolution. Typical applications include:

- Phase identification, standardless quantification, and non-ambient studies on flat powder samples
- Analysis of samples in glass capillaries
- Micro-diffraction

7.1.2 General Setup Procedures (Flat Powder Samples)

Hardware Setup

The line focus of an X-ray tube is traditionally used to collect data for phase measurements of flat specimens because line focus has less divergence than point focus thus allowing a larger area of the sample to be covered.

- After installing an appropriate sample stage, **install either a Programmable Divergence Slit (PDS) or a Fixed Divergence Slit (FDS) optical device as the incident beam PreFIX module.**

When a PDS is used, fixed divergence slit settings of $1/32^\circ$, $1/16^\circ$, $1/8^\circ$, $1/4^\circ$, $1/2^\circ$, 1° , 2° , or 4° are allowed and can be controlled from the *Data Collector* Control Panel. However, the X'Celerator is very sensitive to over irradiation, which effects background characteristics and resolution at low 2θ values. **Consequently, an appropriate incident divergence slit should be selected that permits the sample to intercept the entire spread of the diverging beam at the lowest diffraction angle recorded.**

Additional information: For a fixed goniometer radius, the length L of the sample area irradiated is dependent on 1) the divergence of the X-ray beam and 2) the sin of the angle θ between the incident beam and the sample surface as described by:

$$L = R_o \tan \alpha / \sin \theta$$

where R_o is the goniometer radius in centimeters, and α is the angular aperture of the divergence slit in degrees. A $1/8^\circ$ divergence slit or smaller would be required, for example, to capture the incident beam on a 30 mm sample at $2^\circ 2\theta$ with a goniometer of 240 mm fixed radius (Figure 27).

Note that the use of a PDS in constant illuminated length (variable slit) mode is not currently supported in combination with an X'Celerator in scanning mode in version 1.3x of the *Data Collector* software. Support for variable slits in version 2.0 is unknown at this time.

Specialized Applications

Continued




Figure 27. Divergence slit settings for full illumination as a function of sample length and starting scan angle¹.

¹ For a system with a fixed goniometer radius of 240 mm. For other configurations, see <http://www.geol.uni-erlangen.de/html/xray/vcc.html>

An **Incident Beam Anti-Scatter Slit Holder** (9430 030 86001) can be mounted between the divergence slit and the sample to improve the peak-to-background ratio of the diffractograms (Figure 28).

- The size of this slit should be one step higher than that of the divergence slit, so that when a divergence slit setting of $1/8^\circ$ is chosen, the corresponding anti-scatter slit size should be $1/4^\circ$.




Figure 28. Typical configuration for phase analysis of flat powder samples.

A **Beam Mask** can be inserted after the anti-scatter slit to limit the lateral divergence of the beam (along the goniometer axis) on the sample.

An optional **Beam Knife**, specific for each type of sample stage (Table 3), is recommended for placement just above the sample surface in the center of the goniometer. The beam knife further improves peak-to-background ratios, and is particularly important in scans starting at less than $-6^\circ 2\theta$.

Table 3. Sample stages and corresponding beam knives.

Sample Stage	Beam Knife
9430 030 64601 Reflection/Transmission sample stage	9430 030 64101 Beam knife
9430 030 72601 Stationary Stage for 9430 03018xxx circular sample holders	9430 030 72651 Beam knife for Stage 9430 03018xxx
9430 030 74001 Multi-Purpose Sample Stage (MPSS)	
9430 030 75621 Positioning Stage Z-Tilt-Phi	9430 030 75661 Beam knife for Positioning Stage

The height of the beam knife should be adjusted in such a way that the gap between the knife and the sample is between 0.5 mm and 5 mm. The minimum height of the beam knife should be determined in consideration of the maximum 2θ scan angle as well as the divergence slit setting used (Figure 27). Care should be taken to ensure that the beam is not obstructed by the knife at the highest measurement angles. The exact height positioning of the beam knife is not critical, however, and a setting of ~ 0.5 to 1 mm above the minimum height shown in Figure 2 does not influence the results significantly, especially when using a narrow divergence slit. Detailed instructions for adjustment of the beam knife are provided in the *X'Pert Powder System User's Manual* (Part II, Chapter 12).

Specialized Applications

Continued




Figure 29. Recommended beam knife height for a 240 mm radius goniometer

Soller slits of either 0.02 or 0.04 radians are placed between the source and the divergence slits to collimate the beam in the plane of the sample and to limit asymmetric broadening of the diffraction profile in the low 2θ direction (for angles below $90^\circ 2\theta$).

- A second Soller slit of equal (or greater) size is mounted on the diffracted beam **Anti-Scatter Shield** (9430 030 14001) that mechanically interfaces with the *X'Celerator*.


The Anti-Scatter Shield has a slot positioned after the secondary Soller slit for placement of a **K-Beta filter** or **manual beam attenuator**.

For many powder applications, a Ni filter (9430 031 51031) is used with Cu radiation to cut the intensity of the Cu K-Beta peak to less than 1% of the Cu K-Alpha peak.


- For samples yielding fluorescent radiation (e.g., samples containing high Fe levels when Cu radiation is being used), a graphite Diffracted Beam Monochromator is also recommended.

7.1.3 User Setup and Optics Configuration

1. Open Data Collector enter a user name and password.
2. From the menu, select System Settings > User Management



- Click New



New User

Name:

Full user name:

Access level: ▼

User manager:

Password:


Re-type password:

Buttons: OK, Cancel, Help

- Enter the following information:
 Name: = My Name
 Full user name: = My Full Name
 Access level: = Select "System Owner"
 Check the "User manager" box.
 Password: = password (do not use spaces in your password)
 Re-type Password: = password.
- Click on OK.

Specialized Applications

Continued



- Check that "My Name" is included in the "Name" list and click Close.
The new user is now created. To start actions as this user, select File > Exit to exit from Data Collector and restart Data Collector by double-clicking on the Data Collector icon. Log in with the user name "My Name" and the password "password".
3. Establish communication with the Instrument by selecting '**Instrument**' and then '**Connect**' in the *Data Collector* window.
 4. Select the configuration that specifies which sample stage has been installed, and then press the '**OK**' button to initialize the instrument.
 5. Select the '**Instrument Settings**' tab in *Data Collector* control window and set the generator to 45 kV and 40 mA.
 6. Load samples if an automatic sample stage is installed.
 7. Select the '**Incident Beam Optics**' tab and double click on '**Prefix Module: Progr. Div. Slit Module**' to open the configuration window (Figure 30, right).




Figure 30. Data Collector control window and Incident Beam Optics setup

8. Select the appropriate incident beam module type (e.g., 'Progr. Div. Slit Module' for PDS), and configure each optic by highlighting the corresponding tab and inputting the correct parameter for each installed component. Remember to set the anti-scatter slit to one step higher than that of the divergence slit.
9. Select the 'Diffracted Beam Optics' tab in the Data Collector Control Panel window and double click on 'Prefix Module: X'celerator'. Select the X'celerator module type, and again, configure each optic by inputting the correct parameter for each installed optical component (e.g., Ni filter/beam attenuator, Soller slit, monochromator).
10. Click on the 'Detector' tab, set the detector to 'Scanning mode' in the 'Usage' field, and increase the 'Active Length' to ~3 (Figure 31) . After clicking in any other field to accept this value, the active length will reset to ~2.122, depending on radius of the diffractometer. See 'Modes of Operation' (below) for more information.

Specialized Applications

Continued




Figure 31. Configuration of the detector for 'Scanning mode'.

7.1.4 Modes of Operation

The *X'CELERATOR* is a 'Real-Time Multiple Strip' (RTMS) X-ray detection system that effectively emulates a parallel array of point detectors, each simultaneously collecting data across some finite 2θ range. It can be used in two different modes:

- **Scanning mode** (depicted in Figure 31) is the standard configuration for collecting data over an angular 2θ range. When all parallel detectors are used, the active length is slightly larger than $2^\circ 2\theta$ on a system with a radius of 240 mm. It is possible to reduce the active length of the strips and slightly increase peak resolution by modifying data collection software parameters, but this effect will only be observable when an extremely well-crystalline powder is used. That option is not normally required unless measurement starts at a very low angles (~ 1 degree 2θ). $1/8$, $1/4$ or $1/2$ deg divergence slit settings is recommended depending on the size of the sample and starting angle.
- In **Receiving Slit mode**, the results of all strips are summed. This option is used for texture experiments and for aligning system components and parameters such as X-ray tube height or height of the specimen holder. For optimum sensitivity, we recommend rotating the *X'CELERATOR* 90° on its mechanical interface.

The *X'CELERATOR* can also be used for 2-axis measurements, including the recording of pole figures. In scanning mode, however, only 2-axis measurements where 2θ or $2\theta\Omega$ are one of the scan axes are allowed. Additional information can be found in the *Data Collector* online Help page.




Figure 32. The X'Celerator effectively replaces the detection system of traditional diffractometers by employing more than one hundred thin parallel strips. Just like the traditional theta-theta geometry, the x-ray source and detector are scanned simultaneously. However, during the scan, more reflections are observed simultaneously. Additionally, one particular reflection is followed for a longer time, thus improving sensitivity.

- IMPORTANT!** In addition to configuring the X'Celerator correctly within the Data Collector control panel as described above, you must explicitly specify the mode of operation for the X'Celerator within the 'settings' page of each analysis program employed (Figure 33). When creating or editing an absolute scan program, for example (e.g., in X'Pert Data Collector, 'File > New [or Open] Program'), select 'Settings', double-click on 'Detector: X'Celerator', and then input either 'scanning mode' or 'receiving slit mode' as appropriate. When the mode parameter is set to 'Actual', the X'Celerator will default to the 'Receiving Slit Mode'. When specifying 'Scanning Mode', set the active length to the maximum value of 2.122 as shown in Figure 31. It is best to not leave the Detector specified as "Actual" for "X'Celerator" and "PIXcel". These detectors should be specified explicitly in order for the correct step size and total measurement times to be calculated correctly.



Figure 33. Configuration of analysis program settings for X'Celerator parameters.

Specialized Applications

Continued

7.1.5 Analysis of Samples in Glass Capillaries

The *X'Celerator* can also be used for the analysis of powder samples in glass capillaries using the Capillary Spinner stage. In this situation, the time for recording a useful measurement can be reduced significantly with respect to a proportional detector. When using a momocapillary setup, an additional small anti-scatter shield is typically mounted onto the mechanical interface for optimum performance.

7.1.6 Analysis of Small Samples and/or Micro-Diffraction

For the analysis of small samples or small parts on a larger sample, mono-capillaries are used as the incident beam optic for collimating the X-ray beam to a small, sub-1 mm spot. The *X'Celerator* speeds up the data acquisition considerably and also improves the statistics of the measurement, when compared to a proportional detector.

Examples




Figure 34. Excellent low-angle performance of the *X'Celerator* in the analysis of a Ag-behenate standard. A 1/32° divergence slit was used.




Figure 35. Effect of the beam knife on diffractograms of a poorly crystalline Georgia kaolinite collected using the X'Celerator

7.1.7 PIXcel Detector

The PIXcel is an ultra-fast X-ray detector based on Medipix2 technology, a solid-state technology providing for counting photons with a high spatial resolution (higher than X'Celerator) and a high dynamic range. The PIXcel comprises more than 65,000 pixels, each 55 mm × 55 mm in size. Each pixel has its own counting circuitry, giving rapid readout, a high counting linearity combined with a high dynamic range. The PIXcel I can be used as a linear detector (1D) a point detector (0D) an area detector (2D) or as an imaging detector (3D, only on the Empyrean system). The PIXcel is designed for use with Cu or Co *K_α* radiation.



Specialized Applications

Continued

7.1.8 Applications

The PIXcel can be used in a wide range of diffraction applications. Typical examples of applications are:

- Phase analysis on powder samples.
- Analysis of samples in glass capillaries.
- Texture measurements.
- Parallel beam measurements on thin films and rough samples.
- Micro-diffraction.
- High-resolution rocking curves and reciprocal space mapping.
- Reflectivity.

These applications are all described in detail in the following sections.

7.1.9 Phase Analysis on Flat Powder Samples

Similar to X'Celerator, the PIXcel can be used in 1D scanning mode for powder diffraction applications (phase identification, standardless quantification, crystallography, and non-ambient studies on flat powder samples). The ultra-fast data acquisition makes it possible to optimize on resolution but still have reasonable total measurement times.

A (fixed or programmable) divergence slit with an incident beam anti-scatter slit mounted is used in the incident beam path. An anti-scatter shield is mounted onto the PIXcel.

A beam knife is placed above the sample in order to improve the peak-to-background ratio in the diffractograms.

7.1.10 Analysis of Samples in Glass Capillaries

The PIXcel can also be used in 1D scanning for the analysis of samples in glass capillaries mounted onto a capillary spinner. In this case, the time to record a reasonable measurement can be reduced significantly with respect to a standard proportional detector. The capillary is irradiated by a small beam emitted by a focusing X-ray mirror or a hybrid monochromator. An additional anti-scatter device is mounted onto the PIXcel's anti-scatter shield in order to optimize the data quality.

7.1.11 Texture Measurements

The PIXcel can be used as a point detector (0D or Receiving slit mode) in quantitative texture analysis. The active length of the PIXcel can be set and used as a programmable receiving slit in the Schultz geometry. Alternatively, the PIXcel can be mounted onto a parallel plate collimator for texture measurements in the parallel beam geometry.

7.1.12 Micro-diffraction

When small samples, or small parts of a larger sample need to be analyzed, a mono-capillary is used in the incident beam path to collimate the X-ray beam to a small diameter. Using a PIXcel in the diffracted beam path reduces the measurement time considerably, and improves the counting statistics in comparison to a proportional detector.

7.1.13 Parallel Beam Measurements

In parallel beam applications, such as thin film phase identification, the PIXcel can be used as a point detector (0D, receiving slit mode) mounted onto a parallel plate collimator with the 0D interface for the PIXcel detector.

7.1.14 High-resolution Rocking Curves and Reciprocal Space Mapping

The high dynamic range of the PIXcel makes it a very useful detector for high-resolution applications such as rocking curve analysis and reciprocal space mapping. A manual or programmable beam attenuator is not needed during the measurements. When the PIXcel is used as a point detector (0D, receiving slit mode), it can be mounted onto a rocking curve attachment or onto a triple-axis analyzer; the variable active length can act as a programmable receiving slit.

Note: For Fast Reciprocal Space Map, PIXcel can be used as a line detector (1D scanning or static mode), but the axis of 2D scans must be (2q vs w).

7.1.15 Reflectivity

For reflectivity analysis, the PIXcel can be combined with a parallel plate collimator to be used as a point detector, or it can be used in combination with an anti-scatter slit fixed or programmable) in order to create a beam tunnel (to 0.165 or 0.055 mm). Depending on the intensity of the incident X-ray beam, a beam attenuator may be required, for instance when you are working with an X-ray mirror. For other optical modules, such as a divergence slit or a hybrid monochromator, a beam attenuator is generally not required during the measurements.

7.2 Sample Changer Programs: Data Collector

7.2.1 Sample Changer Programs (for Data Collector)

X'Pert Powder and Empyrean Systems

Hardware Setup

1. Programmable divergence and receiving slits (PDS and PRS) optics are installed on your MPD, and will be controlled through the software.
2. Load individual sample(s) into 15-position magazine(s), and load the magazine(s) into the sample changer. Close the doors.

Create an Analysis Program

3. If the program you need already exists, skip to Step 8. Otherwise, select '**File > New Program > Absolute scan**', to create a new measurement program.

Specialized Applications

Continued




Figure 36. Setup page for analysis (scan) programs.

4. Enter information in the following fields: 'start angle', 'end angle', 'step size', and either 'time per step' or 'total scan time'. All other fields have the correct information by default (e.g., Configuration, Scan Axis, Scan mode).
5. Press the 'Settings' button to configure the analysis program for the desired instrument conditions. Select 'Spinning', and set the speed to approximately match the analysis step time (e.g., set the spin rate at 1 revolution sec^{-1} if the stepping rate is 1 sec step⁻¹.)
6. Select 'Incident Beam', and set the 'Prog. Div. Slit' to 'Fixed', and 'Angle' to '1°'. On the 'Diffracted Beam Optics' tab, select 'Prog. AS Slit', set 'Usage' to 'Fixed', and 'Angle' to '1°'. Next to 'Receiving slit', select 'Prog. Rec. Slit' and set the 'Height' to '0.3 mm'. Other settings, including automatic irradiated length, can also be chosen as needed.
7. From the command line menu select 'File > Save as' and name the program. It might be helpful to name it according to the scan parameters, e.g., 25-75deg pt02 step 1sec-step, so that the scan can be easily reused for different analyses.

Create a Sample Changer Batch Program

8. Select 'File > New Program > program type: Sample (changer) batch'. Press the 'Insert' button and at the top of the new window, select 'Browse' to point to the desired analysis program. Select this program and enter a suitable 'Dataset name' – this will be the file name and it must be unique. Finally, at the bottom of the window, you must specify the sample position in the magazine (i.e., number 1-15).
9. Continue to insert the program(s) to be run for each sample in the changer. Press the 'Sort' button to organize the programs in order of sample changer positions, and press 'Check' to ensure that there are no mistakes with the file names. The total time for the batch is also calculated and shown.





Figure 37. Setup page for Sample Changer (batch) programs.

10. Save the Sample (changer) batch by going to selecting 'File > Save as' and naming the program.

Run the Sample Changer Batch Program

11. Click on 'Measure > Program > Sample (changer) batch' and select the program used in Step 10. Click 'OK' to start the analysis



Specialized Applications

Continued

7.3 Non Ambient Scans

7.3.1 Collecting Non-Ambient Scans Using *Data Collector*

High-Temperature Experiments with the Paar HTK-1200N

Hardware Setup

- Dismount the current sample stage. If there is a sample changer installed, ensure that it is “parked” on the locking ring.
- Carefully mount the HTK-1200N. Initially, bring the camera in while tilting it back at an angle toward the top of the stage interface to hook under the support lever. Make sure the camera rests on the alignment blocks.
- Tighten all bolts, being careful not to over torque them.
- Mount the diffracted beam optics of choice, depending on the material to be analyzed – use either the mirror/PRS with Xe miniprop detector, or the *X'Celevator* detector.
- Remove the beam mask from the incident optic. The HTK-1200N is self-masking by virtue of the width of the windows.

Sample Prep

- Open the HTK-1200N housing by loosening the four black knobs on bottom until the sample stage drops down out of the camera.
- Unscrew the sample holder from post carefully.
- The sample should be as fine-grained as possible (i.e., the Paar manual suggests 10 μm). An internal standard such as gold powder should be used for absolute temperature determination based on lattice parameter.
- Reattach the holder to the post and carefully insert the sample stage into the camera, taking care to ensure that it rides up along the two metal guide posts. Tighten the black knobs to lock in the sample stage in place. If HTK-1200N has a spinner stage, activate the spinner switch on the cable attached to the camera.

Preparing the Experiment

- Open the *Data Collector* program.
- Click on '**Instrument > Connect**' and select the configuration that defines the HTK-1200N camera. Press the '**OK**' button, and press '**OK**' when the advice list comes up. You will now see the Control Panel window.
- Select the '**Incident Beam Optics**' tab and change the divergence slit to the desired setting. Note that 1° is a 'standard' fixed slit setting; see Section 7.1 for more information).
- Select the '**Diffracted Beam Optics**' tab and set the diffracted beam optics as follows: for mirror/PRS, set the PRS according to the resolution required (0.3 mm is a medium resolution setting, 0.1 is a high resolution setting that will dramatically decrease intensity). The *X'Celevator* detector should be configured for '**Scanning mode**' with an active length of ~3° (see Section 7.1 for more information).
- Select the '**Instrument Settings**' tab and set the generator to 45 kV and 40 mA.
- Press '**OK**'.

Installing the Vacuum System

- Vacuum is necessary for all experiments carried out above 300°C!! Remove the vacuum feed-through covers, and lightly grease the O-ring if needed before connecting.
- When starting operation of the HTK-1200N:
 - Start vacuum system
 - Temperature control unit – press the '**Heater on**' button, which will also start water flow cycle

Determining Proper Measurement Conditions

1. Determine the optimum measurement conditions for a sample, using the procedures outlined in the previous sections.
2. Once the desired measurement parameters (e.g., step size, scan range, etc.) have been determined, create (or modify) a new scan program as outlined in Section 4.1.1 and save the program.

3. Create a batch program for non-ambient analysis by selecting 'File > New Program > Non-ambient program'. Insert the following information into the non-ambient program:
 - an *initial temperature* setting,
 - a *dwel time* at that temperature,
 - a *measurement (absolute scan) program*.
 Keep inserting these three parameters in that order to build up your profile (Figure 38).
4. Press the 'Profile' button to display a schematic of the non-ambient program steps that have been input into the program (Figure 39).
5. Select 'File > Save as' and enter a name for the analysis program.




Figure 38. Creating a non-ambient analysis program.




Figure 39. Graphic display of the non-ambient program steps.

Specialized Applications

Continued

Measuring the Sample

1. Select '**Measure > Program > Non-ambient**' and select the program name used in Step 5 (above).
2. Enter a '**Data Set Name**' and '**Description**'.
3. Press the '**Start**' button.

Completing the Experiment

VERY IMPORTANT!!

- After the analysis and the camera has cooled to ambient temperature (~25°C), you must vent the vacuum **before** turning off the vacuum pump. **This prevents backstreaming of pump oil into the chamber!!** Vent slowly to prevent disturbing the sample.

7.3.2 High-Temperature Experiments with the Paar XRK-900

Hardware Setup

1. Make sure sample changer is "parked" on its locking ring. Dismount current stage.
2. Insert the set of long bolts a few turns into the outer set of screw on the stage interface. Make sure the washers are close to the bolt head before mounting the XRK camera.
3. Carefully mount the XRK by tilting it back and by coming in at an angle toward the top of the stage interface to hook under the support lever. Then slip the large holes on the XRK over the bolts you installed, making sure the camera rests on the alignment blocks.
4. Working on one side of the camera at a time, place and hold the boomerang shaped metal piece under the bolt head and washers.
5. Holding this piece in place with one hand, use the other hand to tighten the bolts with the 10mm wrench provided. Repeat this process on the other set of bolts, finish tightening all bolts, being careful not to over torque the bolts.
6. Remove beam mask from incident optic – the XRK is self-masking by virtue of the width of the windows.

Sample Prep

7. The sample stage is held inside the camera by 4 knurled knobs under the stage. Loosen these knobs to lower and remove the stage and to load a sample. Reinsert the stage and finger-tighten the knobs.
8. Follow the XRK instructions for attaching a reaction gas line and modulating the flow through the sample.

Software Setup

- Open the *Data Collector* program.
- Click on '**Instrument > Connect**' and select the XRK reaction chamber configuration. Press the '**OK**' button, and then again press '**OK**' when the advice list comes up. You will now see the Control Panel window.
- Select the '**Instrument Settings**' tab and double-click on the '**Tension**' setting. Set the generator to the desired power level.
- Select the '**Incident Beam Optics**' tab and change the divergence slit to the desired setting. Note that 1° is a 'standard' fixed slit setting (see Section 7.1 for more information).
- Select the '**Diffraction Beam Optics**' tab and set the diffracted beam optics as follows: for mirror/PRS, set the PRS according to the resolution required (0.3 mm is a medium resolution setting, 0.1 is a high resolution setting that will dramatically decrease intensity). The *X'celerator* detector should be configured for 'Scanning mode' with an active length of ~3° (see Section 7.1 for more information).

Installing the Vacuum System

- Vacuum is necessary for all experiments carried out above 300°C!! Remove the vacuum feed-through covers, and lightly grease the O-ring if needed before connecting.
- Maintain following sequence when starting operation of the XRK:
 - Vacuum pumps
 - Vacuum meters
 - Cooling water cycle
 - Temperature control unit – press the 'Heater on' button, which will also start water flow cycle

Determining Proper Measurement Conditions

1. To determine the optimum measurement conditions for a sample, collect a Manual scan using the procedures outlined in Section 4.1.1.
2. Once the desired measurement parameters (e.g., step size, dwell time, scan range, etc.) have been determined, create (or modify) a new scan program as outlined in Section 4.1.1 and save the program.
3. Create a batch program for non-ambient analysis by selecting 'File > New Program > Non-ambient program'. Insert into this non-ambient program:
 - an *initial temperature* setting,
 - a *dwell time* at that temperature,
 - a *measurement (absolute scan) program*.Keep inserting these three parameters in that order to build up your profile (Figure 38).
4. Press the 'Profile' button to display a schematic of the non-ambient program steps that have been input into the program (Figure 39).
5. Select 'File > Save as' and enter a name for the analysis program.

Measuring the Sample

6. Select 'Measure > Program > Non-ambient' and select the program name used in Step 5 (above).
7. Enter a 'Data Set Name' and 'Description'.
8. Press the Start button.

Completing the Experiment

VERY IMPORTANT!!

- After the analysis and the camera has cooled to ambient temperature (~25°C), you need to vent the vacuum **before** turning off the vacuum pump. **This prevents backstreaming of pump oil into the camera!!** Vent slowly to prevent disturbing the sample.

7.4 CHC Humidity Chamber Measurements

7.4.1 Humidity Chamber Experiments

Hardware Set up

1. Install the Anton Paar CHC stage



Specialized Applications

Continued

NOTE: Attach the sign printed on page 4-20 to the XRD to protect all users from hot components!


2. Prepare a powder sample in the standard sample holder that comes with the stage and insert the holder into the stage.



3. Close the door of chamber and lock it down.
4. Power on the Neslab/Jalabo water bath, the MHG humidity generator, and the TCU-50 temperature controller.
5. Set the above three devices according to the experiment chart shown below. To adjust the setting on the Neslab/Jalabo water, press and hold the display button, and then use the course and fine knobs to set the desired temperature.

Experiment Design

Open the MHG Method Editor Program




MHG Method Editor Program

- portional Detector Xe




Sample Stage:


Five-axes cradle (chi,ph)



First select the segment type to put in a ramp of the humidity from 0 to 80% RH and put in a time of 10 minutes to get to that RH and hit "add segment."



Next add another segment of step for 10 minutes which will keep the RH at 80% for 10 minutes.



Now ramp of the humidity from 80 to 0% RH and put in a time of 10 minutes to get to that RH and hit "add segment."

Specialized Applications

Continued



Then click on the save method button and save the method you have just created.

User Setup – XRD software

Measurement Programs

To control the MHG temperature and humidity controller it is necessary to prepare some measurement programs in the PANalytical Data Collector in the way as described below.

- The Stationary Measurement program is used as a dummy measurement to transfer the desired non-ambient settings via the scripting file SETMHG.exe to the MHG temperature and humidity controller.
- The Absolute Scan program specifies your result measurement used with your sample under investigation.
- The General Batch program specifies the sequence of non-ambient settings necessary for the overall result measurements.

Stationary Measurement

1. Select menu item *File, New Program* and select program type *Stationary measurement*.
2. Ensure that in the drop-down listbox the configuration *AP CHC Plus + MHG Humidity Generator* is selected.
3. Change the time to 1 second, see Figure 6.
4. Save this program under the **pre-defined name** *SetMHG*.

NOTE: Do not use any other names because the delivered script works only correct with this naming convention.




Figure 40. Stationary Measurement program for a dummy program with its prescribed parameters for the use in the General Batch program.

Absolute Program as Result Measurement

1. Select menu item *File, New Program* and select program type *Absolute Scan*.
2. Ensure that in the drop-down listbox the configuration *AP CHC Plus + MHG Humidity Generator* is selected.
3. Ensure that for line detectors of type X'celerator or PIXcel in the *Settings* window the appropriate detector mode is applied, e.g. *Scanning line detector* mode.
4. Change the default measurement parameters of scan axis, start angle, end angle, step size, and time per step to the desired ones, see for an example Figure 41.

Save this result measurement program under your desired name appropriate for your sample, for example *MHG Result Humidity* as shown in Figure 41.




Figure 41. Example for a result measurement program of type *Absolute* for line detectors.

General Batch Program

The last measurement program combines all the necessary dummy programs simulating the non-ambient settings with the desired result measurements, see Figure 46.

1. Select menu item *File, New Program* and select program type *General Batch*.
2. Ensure that in the drop-down listbox the configuration *AP CHC Plus + MHG Humidity Generator* is selected.
3. Select button *FILE NAME SETTINGS* and **switch off** the options *Use the temperature in the file name* and *Use sample ID in the file name*, see Figure 42.
4. Create the desired batch program by selecting button *INSERT MEASUREMENT PROGRAM* and selecting the stationary program *SetMHG*, see Figure 43.
5. Enter in field *File name suffix* the desired values for the non-ambient settings with the following necessary convention: "*<T-value> <underscore> <humidity value> <underscore> <flow value>*", e.g. *40_50_300*, see Figure 43.
This means that the result temperature should be 40 °C and the result humidity is 50 %RH with a flow rate of 300 ml/min.

Specialized Applications

Continued

NOTE: do not change the field *File name prefix* because the name of the stationary program *SetMHG* must be preserved, see above.

6. Select button *Insert Timer Settings* and enter a wait step with at least 5 min to ensure that the correspondent temperature and/or humidity step is properly executed before starting the result measurement, see for example Figure 44.
NOTE: a better value can be derived by earlier made test runs to check when the time is arrived the temperature and the humidity will be stable inside the non-ambient chamber.
7. Select then again button *INSERT MEASUREMENT PROGRAM* and select the absolute scan program used for the result measurement, e.g. *MHG Result Humidity*, see Figure 45.
8. Enter in field *File name suffix* the specified values for the non-ambient settings using the same convention as above "*<T-value> <underscore> <humidity value> <underscore> <flow value>*", e.g. *40_50_300*, see Figure 45.
In this way the name of the xrdml result measurement contains the correspondent non-ambient settings values for proper reporting.
9. Repeat steps 4 to 8 and build up your non-ambient program, see for example Figure 46.

NOTE: When preparing the General Batch program take into account the temperature and humidity limits.

For the Anton Paar CHC Plus non-ambient stage with MHG humidity generator these are:

0.0 £ T (°C) £ 80.0 and 0 £ Hum (%RH) £ 100.




Figure 42. Batch File Name Settings window with indicated switched off functions.




Figure 43. Measurement Program Step window with necessary specified non-ambient settings.




Figure 44. Timer Settings window with minimum wait step.




Figure 45. Measurement Program Step window with necessary specified settings for a result measurement and indicated non-ambient settings for reporting in xrdml file name.

Specialized Applications

Continued




Figure 46. Example of a General batch program simulating the non-ambient parameters for the sample stage AP CHC Plus + MHG Humidity Generator.

Optics/Instrument Setup

1. Select the Goniometer & Sample Stage tab.
2. Press the button 'Change kV-mA' to set the generator to 45 kV and 40 mA.
3. Press 'OK'
4. The optics can either be set inside your absolute scan program (See page 27 "Prepare Absolute Scan" of this document), or by doing the bookkeeping in the Control window. To set the optics from the Control window, select the Incident beam tab (See page 25).
5. Press the 'Change Optics' button and set the divergence slit to 1° fixed. For low-angle work (<5°), select ½ deg. divergence or smaller for fixed divergence slit depending on sample length and lowest starting angle. For systems with PDS select Automatic mode with a 10 mm length or whatever matches the sample length.
6. Select the 'Diffracted Beam Optics' tab and set the diffracted beam optics as follows: set FRS or PRS according to the resolution required (0.3 mm is a medium resolution setting, 0.1 is a high resolution setting that will dramatically decrease intensity). The X'Celerator detector should be configured for 'Scanning mode' with an active length of ~3° (see Section 7.1 for more information).

Ending the Experiment

It is very important to follow the subsequent instructions at the conclusion of the experiment, particularly because when the sample stage is heated, possible burn injury may result if the equipment is mishandled!

1. Reset your equipment as follows:
 - TCU50 – set sample stage to ambient (20° – 25°C)
 - Set the Neslab/Jalabo to 25°C
 - Set MHG humidity generator to 35°C. Keeping the MHG on at all times at ~ambient temperature will prevent dry-out and the need for a ~24 hour stabilization period after re-start.
2. Allow these units to reach the above temperatures, particularly the sample stage. Once they have reached these ambient temperatures, it is safe to remove the sample from the chamber.

7.5 Stress

Stress Measurement Acceptance Testing using the Data Collector

System Configuration

Incident beam optics:

- PreFIX crossed slits collimator (knob adjusted).

Diffacted beam optics:

- Programmable Receiving Slit (PRS)
- Diffracted Beam Mask interface for texture.
- Programmable Anti Scatter Slit (PASS) or Fixed Anti Scatter Slit (FASS) can be mounted to the PRS.

Detectors:

- Proportional Detector Xe.

Sample Stage:

- Five-axes cradle (chi, phi, x, y, z) for 240 mm radius

To carry out this test using the Data Collector, a configuration must be available or created in which the X-ray lens and/or the Mono capillary and or the Crossed Slits Collimator are available.


1. Open the Data Collector program.
2. Select or open the appropriate configuration and go on-line; the Control window pops-up.
3. Check if all optical items are selected for the measurement and that the settings are correct.

Specialized Applications


Continued

Modified Stress Setup


4. Set Chi = 45 deg.




5. Select the User settings pull-down menu and then select the Sample Offsets menu. Create a sample offset of 45 ° for the chi-axis and press OK.



- a) With mouse point over "Chi" axis press right mouse button and choose "Change". Change Chi axis = "45" to "0". Sample offset for "Chi" will change from "0" to "-45" deg.





- b) Note that when Okay button is pressed the Chi axis is recalibrated to "0".



Specialized Applications

Continued


- From Pull-down menu Select, File, Create a new program and select Stress measurement.




- a. Enter the following scan parameters:

Item	Parameter
Tilt Axis	Chi
Tilt range	Positive + negative
Maximum Tilt	45
No. of Tilt steps/Phi	2
Scan Axis	gonio
Scan mode	continuous
Range	4°
Step size	0.04°
Time per step	10.0 seconds


b. Changes shown below.



8. Press the Scan table key
 - a) Press Yes button to view the Scan table.



b) Scan Table.



No.	Psi (°)	Sin2Psi	Phi (°)	2Theta (°) Start	2Theta (°) End	2Theta (°) Step Size	No. of steps	Time per step (s)	Net time per step	Pre-set counts	Scan time (h:m:s)
1	45.00	0.5000	0.0	142.0200	145.9800	0.0400	99	10.00	10.00	10000	00:16:32
2	0.00	0.0000	0.0	142.0200	145.9800	0.0400	99	10.00	10.00	10000	00:16:32
3	45.00	0.5000	0.0	142.0200	145.9800	0.0400	99	10.00	10.00	10000	00:16:32

9. Check if the correct start and end angles are filled in.
10. Press the OK key.
11. Save the program, for example, as: **Stress acceptance test**.
12. Execute the **Stress acceptance test** program.
13. After completion open the Data Viewer program and select the scans (combined graph)
14. Select the tools pull-down menu and put a check mark in front of Spectral lines.
15. Determine the peak positions by shifting the line and check if the three peaks are overlapping each other within **0.01°**. If correct then this completes the chi-stress measurement for the X-ray lens and /or mono-capillary.

Specialized Applications

Continued

7.6 Texture

Texture Measurement Acceptance Testing Texture Measurement with Crossed Slits Collimator

System Configuration

Incident beam optics:

- PreFIX crossed slits collimator (knob adjusted).

Diffracted beam optics:

- Programmable Receiving Slit (PRS)
- Diffracted Beam Mask interface for texture.
- Programmable Anti Scatter Slit (PASS) or Fixed Anti Scatter Slit (FASS) can be mounted to the PRS.

Detectors:

- Proportional Detector Xe.

Sample Stage:

- Five-axes cradle (chi,phi,x,y,z) for 240 mm radius
1. Install the Crossed Slits Collimator.
 2. Install the PRS and remove the Soller slits. If available, install the beam mask holder of the texture kit and insert a beam mask 10.
 3. Mount the textured Cu sample in the center of the Five-axes cradle (chi,phi,x,y,z) for 240 mm radius. If strations are visible on the Cu sample surface, the strations correspond to "Rolling Direction". Mount the Cu sample so that the strations are perpendicular to the direction of incident beam
 4. Insert the Ni filter (Fe for Co) into the filter holder of the PRS (or into the Crossed Slit Collimator).
 5. Set the slits of the Crossed Slit Collimator to **1 × 1 mm**.
 6. Select **40 kV** and **50 mA** on the generator for Cu. (or 35 kV and 50 mA for Co)
 7. Leave TDS and enter the Data Collector program.
 8. Select or open the appropriate configuration and go on-line; the Control window pops-up.
 9. Check if all optical items are selected for this measurement and that the settings are correct.
 10. Enter the **File** from the main menu and select from this menu the **New Program** option.
 11. Select from the selection menu **Program type: Texture measurements**.
 12. Press the **OK** button; the **Prepare Texture measurement** input menu is shown.
 13. Enter the following program parameters:
 - Configuration:** Select the appropriate configuration for the Crossed Slits Collimator configuration. Also select the required diffracted beam path. (from the list, if programmed in the configuration)
 - Measurement type:** Single pole figure measurement.
 - 2Theta:** 74.15° (Cu) or 88.865° (Co)
 - hkl for ODF:** 2 2 0

14. Click on **Settings**; the Settings input menu is opened.
 - a. Double click on Sample Stage and set the Oscillation to 5mm.
 - b. **"Incident beam path"**:
 - PreFIX module: actual.
 - "Anti Scatter Slit" : none.
 - "Beam Attenuator": none."Divergence Slit": Crossed Slit Collimator, 1 mm.
 - "Filter": none. If required Ni filter (Cu) or Fe (Co)
 - "Capillary Optic": none.
 - "Mask": Crossed Slit Collimator, 1 mm.
 - "Mirror": none.
 - "Monochromator": none.
 - "Soller Slit": none.
 - c. **"Diffracted beam path"**:
 - "PreFIX module: actual.
 - "Anti-scatter slit": Prog. AS Slit, Fixed, Angle: 4° (only if present, else to "none")
 - "Beam attenuator": none.
 - "Collimator": none.
 - "Detector": PW3011/20 [Miniprop. large window]
 - "Filter": Ni or Co if not used in the incident beam path.
 - "Mask": Mask Fixed 10.
 - "Mirror": none.
 - "Monochromator": none.
 - "Receiving slit": Prog. Rec. Slit, Height = 2 mm.
 - "Soller slit": none.
15. Click the **Apply** button to activate the selected settings.
16. Click the **OK** button to return to the **Prepare Texture measurements** menu.
17. Click the **Edit** button to activate the Prepare Pole figure.
18. Enter the following parameters:

Item	Parameter
Scan mode	Continuous
Phi	Start angle = 0°, End angle = 360°, Step size = 3°, Time per step = 0.65
Chi	Start angle = 0°, End angle = 84°, Step size = 3°

19. Click the **OK** button to return to the **Prepare Texture measurements** menu.
20. Select the **File** from the **Main** menu and then the **Save As** option; the **Save Program As** input menu appears.
21. Enter in the **"Name"** field: "Texture" and in the **"Description"** field: "Texture measurement with Crossed Slits Collimator and Cu sample".
22. Press the **OK** button to save the program.

7.6.1 Starting the Texture Measurement

1. Select **Measure** from the main menu and then **Program**; the Open Program selection menu pops-up.
2. Select from this menu the **Texture measurements** and from the list the **previous stored program**.
3. Enter in the Start menu the **Project name** and the **Data set name**.
4. Check if all positions are still as were programmed or set in the control window.
5. Press OK to start the Texture measurement.
If an earlier measurement was already stored with the same name then a message box will appear. Select the appropriate button to start the measurement.
6. After completion check if the system moves to its rest positions, if filled in.

Specialized Applications

Continued




Figure 47. Typical Pole Figure Plot using the X'Pert Texture Software

7.7 Thin Film Grazing Incidence Analysis

7.7.1 PreFIX Systems: Optics and Sample Changers

Thin Film / Grazing Incidence Analysis

Hardware Setup

The line source of the X-ray tube is used for thin film analysis, either with a PDS as the incident optic, or preferably with a parabolic X-ray mirror to convert the divergent X-ray beam to an intense, monochromatic, quasi-parallel beam (Fig. 5.9). A $1/32^\circ$ slit must be inserted in the mirror. An appropriate beam mask and an attenuator foil must also be inserted into the mirror, depending on sample characteristics.

A parallel plate collimator must be mounted on the diffracted beam side, with a 0.04 rad Soller slit inserted.




Figure 48. Comparison of a) focusing geometry, b) parallel beam geometry with a PDS incident optic, and c) parallel beam with an incident mirror.

User Setup

1. Open *X'Pert Organizer* and select an appropriate Users and/or Projects.
2. Select **'Modules > X'Pert Data Collector'** or the *Data Collector* button on the tool bar.
3. Select **'Instrument > Connect'** and input the correct configuration based on the stage being used.

Optics Setup

4. Select the **'Incident beam optics'** tab and double click on any item in the list. Change the optics to reflect the current instrument setup. The mirror will be selected as the *PreFix* module. Input the correct values for the Soller slits, mask, and other components.
5. Select the **'Diffracted beam optics'** tab and input **'Parallel plate collimator'** as the *PreFix* module. The Soller slits will be set to 0.04 rads, the detector will be mini prop with the wavelength set to K-alpha
6. Select the **'Instrument'** tab and set the generator to 45 kV and 40 mA.

Measurement Programs

7. Create a analysis program by selecting **'File > New Program > Absolute scan'** from the command line.
8. Enter the following information:
 - 'Scan axis': **'2theta'**
 - 'Other gonio angle': set to omega, and select a grazing angle between 0.5 and 3 degrees, depending on how shallow or deep you wish to penetrate into the sample.
 - 'Scan mode': **'Continuous'**

Specialized Applications

Continued

- 'Start and end angle(s)': Enter values based on where the reflections are expected for the material being studied.
- 'Step size' and time per step should be appropriate for the material you are analyzing (see manual scan instructions for how to determine these)




Figure 49. Input parameter page for a grazing incidence experiment.

9. Select 'File > Save as' and enter a name and description for this scan program and then push 'OK'.

Collecting Measurements

10. Select 'Measure > Program > Absolute scan' and input the program name used in Step 2.
11. Enter a 'Data set name' and 'Description'. Press 'Start'.

7.8 Reflectivity

7.8.1 Reflectivity measurement using standard reference thin film sample

System Configuration

The System must be configured as shown in the figure below.




Figure 50. Basic Thin Film Configuration with Parallel Plate Collimator (PPC)

Incident Beam Optics

- Preferably a Programmable Divergence Slit (PDS) with a Programmable Beam Attenuator Ni 0.125 mm (PBA) (for Cu radiation) mounted
or
- A Fixed Divergence Slit (FDS) with Programmable Beam Attenuator Ni 0.125 mm (PBA) (for Cu radiation) mounted.
- Use a set of Low Angle Slits in combination with a Fixed Divergence Slit.

Diffracted Beam Optics

Three types of Parallel Plate Collimators can be used, they are:

- 0.09 deg.
- 0.18 deg.
- 0.27 deg.

Detector

- Proportional Detector Xe, PIXcel3D Detector Single or X'Celerator Detector.

When using a PIXcel3D detector or X'Celerator a 0-d interface must be fitted to the detector.

Sample Stage

- Stage for flat samples. (supplied with the system)

For other sample stages that can be used refer to the relevant chapters in the Empyrean User's Guide.

7.8.2 Basic Thin Film Configuration Acceptance Test

The basic thin film configuration acceptance test can be carried out with or without a Flat crystal monochromator fitted. The purpose of this test is to measure a "fringe pattern" totally reflected from the Thin Film sample (Cr).

If an attenuator is fitted, it will be automatically switched out of the beam path during the test.

Related information


Specialized Applications

Continued

7.8.3 Testing a Thin Film Configuration with a Proportional Detector

A Flat diffracted beam monochromator is optional for better background suppression if fluorescence from the sample is expected.


1. Install the Parallel Plate Collimator with detector in the diffracted beam PreFIX holder and insert the **0.1 reflectometry slit**.
2. Select the appropriate diffracted beam optics for the parallel plate collimator.
3. Set the HT generator to **45 kV** and **40 mA**.
4. Set the Divergence to $1/32^\circ$ or insert $1/32^\circ$ divergence slit.
5. Check whether the **0.04 rad.** soller slit and the **10 mm** beam mask are installed.
6. Insert a **Ni filter** into the filter holder if no Flat Crystal monochromator fitted.
7. Place the **Cr Thin Film reference sample** into the Stage for flat samples.
8. Open the X'Pert Datacollector.
9. Select or create a configuration that has the Parallel Plate Collimator(s) listed.
10. Carry out a manual scan using a Range of **0.7 degrees**, a Step size of **0.005 degrees** and a Time per step of **0.2 seconds with the system at a 2theta angle at 1.25 degrees**.
11. After completion the fringe pattern as shown in the figure below should be visible.



12. Repeat this procedure in case other parallel Plate Collimators are part of the configuration. The intensity will be lower but the fringe positions must be at the same location within **0.01° 2theta**.

Testing a Thin Film Configuration with a PIXcel^{3D} Detector


1. Install the Parallel Plate Collimator with PIXcel^{3D} detector (using 0D interface) in the diffracted beam PreFIX holder and insert the **0.1 reflectometry slit**.
2. Repeat measurement as stated above in Proportional detector section.
3. After completion a fringe pattern, as typically shown in the figure, is plotted.



4. Repeat this procedure in case other parallel Plate Collimators are part of the configuration. The intensity will be lower but the fringe positions must be at the same location within 0.01° 2θ .

Testing a Thin Film Configuration with a X'Celerator Detector

1. Repeat measurement as stated above in Proportional detector section.
2. After completion a fringe pattern, as typically shown in the figure, is plotted.



3. This completes the functional test of the Parallel Plate Collimator(s).


Specialized Applications

Continued

7.9 Data Base Conversion and Installation

7.9.1 Installing and Converting Current ICDD Databases (PDF-2, PDF-4+, PDF-4/Organics)

1. Insert the CD and run the Setup executable (.exe).
2. Once the Installation is complete, you will need to obtain a license by running the **Registration Wizard**.
 - a. Click **Start**.
 - b. **Go to All Programs**
 - c. Find the folder for the product you installed
 - d. Choose Register ICDD Product.
 - e. Select: Obtain a registration key for your product from ICDD.
 - f. Fill out required information and send the information whichever way you prefer.
 - g. Once you receive the Registration Key, re-run the Registration Wizard.
 - h. Select: Enter a registration key obtained from ICDD.
 - i. Enter the required information.
3. Open **HighScore Plus**
4. Go to *Customize > Manage Databases*
5. If the database you have installed does not appear at the top under the column *Database name*, click **Search for Databases** (see below).



6. Once the database appears, click the column use, and then press the *Convert* button (see below)



7. Once the database is converted (can take up to a few hours) the same database will appear, but will now be a HighScore (Plus) V3.X database. You will want to check Use for this database, and un-check Use for the old, Sybase database which you just converted.

Note: ICDD databases can be used directly without a conversion. It is **recommended** to convert these databases for a much faster access by the HighScore software.

7.9.2 Installing the PAN-ICSD

The PAN-ICSD database contains inorganic crystal structures and is derived from the FIZ ICSD database. It works only in combination with HighScore (Plus) and completes PDF-2 patterns with structure information.

1. Insert the CD and run the Setup executable (Setup.exe).
2. The presence of a PAN-ICSD database is automatically detected and it shows up in the database manager.
 - a. Check the box Use to use this database.

7.9.3 Converting Legacy Databases (PDF-2, PDF-4+, PDF-4/Organics)

These are old reference database formats used by former Philips or PANalytical software.

The following databases can be converted for use in HighScore:

APD1700 format

PC-APD format

PC-Identify/X'Pert Graphics and Identify format


PDF1 reference database in PHILLIPS format (NOT the original PDF-1)

PDF-2 ICDD Products


Note: "Converted" databases in the previous version 2.0 PANalytical database format:

Just keep on using these converted databases in version 3.0 of HighScore. There is no need nor benefit in converting these databases into the version 3.0 database format.

To Convert these legacy databases, you will need to run the external program, **RefDBConversion.EXE**, which can be found in the HighScore installation folder.



Double click on this icon, and a new window will appear



Choose the format of the database you wish to convert, and then enter the remaining required information and click **Start Conversion**.